

PAWTAPE

Synthetic Veterinary
Adhesive Tape

Milliken



PawTape is a non-stretch, hand-tear synthetic adhesive tape featuring fun paw prints. It is the first adhesive tape designed for small animal use to be engineered with a non-stretch, moisture-stable substrate.

Product Applications

- Immobilize wounded areas
- Splinting, strapping, and patient positioning
- Hold other dressings in place
- Protection and support
- Synthetic alternative to cotton tape



How to Tear



Pinch, Pinch, Rip!

Pro Tip:

PawTape is **STRONG** - tear quickly for cleanest results!



How to Apply



After cleaning and bandaging the area of application, place the adhesive tape half on the bandage and half on the skin.



Apply as necessary.



Hand-Tear or cut with scissors and place hand over tape for upwards of 30 seconds to slightly warm the adhesive.

Features:

Strong

33 lbs. tensile strength

Synthetic

cotton alternative

Easy Hand Tear

no scissors needed

Consistent support

will not constrict

Water-resistant

Ordering Information:

Catalog #	Size	Rolls/Box	Boxes/Case	Rolls/Case
291-005-096N	½" x 10 yds.	24	4	96
291-010-048N	1" x 10 yds.	12	4	48
291-020-024N	2" x 10 yds.	6	4	24