

NUSTAT

Topical Hemostatic Dressing



NUSTAT Flex dressings help address patient blood loss through hemostatic fiber technology. The specialty fibers within the product result in a flexible dressing that conforms to the patient. NUSTAT Flex dressings are applied externally with mechanical compression to temporarily control bleeding.



INDICATIONS FOR USE

NuStat is a single-use hemostatic wound dressing applied externally with mechanical compression to temporarily control bleeding in:

- Lacerations
- Punctures
- Abrasions
- Surgical wounds (operative, postoperative, dermatological, etc.)
- Traumatic injuries



BENEFITS

- Fast - Quickly gain control of bleeding
- Flexible - Conforms to target wound area
- Durable - Will not tear when wrapped
- Cost Effective



Milliken™

MILLIKEN HEALTHCARE PRODUCTS, LLC

ORDERING INFORMATION

Size	Code	Packaging
3" x 48"	3000097317	10 ea/cs
4" x 48"	3000097320	10 ea/cs
12" x 12"	3000097333	10 ea/cs
6" x 60"	3000097331	10 ea/ca

Size	Code	Packaging
Trauma Pad 8" x 12"	3000097832	10 ea/cs

Two Easy Ways to Order:

1. Preferred Distributor
2. Direct from Milliken

medproducts@milliken.com | healthcare.milliken.com | 1-866-491-6556